

## Même®

textured, silicone gel-filled  
mammary implants, round shape,  
available in four projections

# POLYTECH

Health & Aesthetics











The silicone gel-filled Même® implants with textured surface have a round base and a central projection.





The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections: Low Profile (TLP), Moderate Profile (TMP), High Profile (THP) and Extra High Profile (TXP).





The new curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the infra-mammary fold.

# Même® , CE 0483

Même® TLP – Low Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
	A	C	B	D	
20624-100*	93	93	25	54	<b>100</b>
20624-115	98	98	26	57	<b>115</b>
20624-135*	103	103	27	60	<b>135</b>
20624-150	108	108	28	62	<b>150</b>
20624-170*	113	113	28	65	<b>170</b>
20624-195	118	118	29	68	<b>195</b>
20624-215*	123	123	30	71	<b>215</b>
20624-245	128	128	31	73	<b>245</b>
20624-270*	133	133	32	76	<b>270</b>
20624-300	138	138	33	79	<b>300</b>
20624-330*	143	143	34	82	<b>330</b>

Même® TMP – Moderate Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
	A	C	B	D	
20625-110	88	88	31	56	<b>110</b>
20625-130*	93	93	32	59	<b>130</b>
20625-150	98	98	34	62	<b>150</b>
20625-175*	103	103	35	65	<b>175</b>
20625-200	108	108	36	68	<b>200</b>
20625-230*	113	113	37	71	<b>230</b>
20625-260	118	118	39	74	<b>260</b>
20625-290*	123	123	40	77	<b>290</b>
20625-325	128	128	41	80	<b>325</b>
20625-365*	133	133	42	83	<b>365</b>
20625-405	138	138	44	86	<b>405</b>

Même® THP – High Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
	A	C	B	D	
20626-140	88	88	38	61	<b>140</b>
20626-165*	93	93	39	64	<b>165</b>
20626-190	98	98	41	67	<b>190</b>
20626-215*	103	103	43	71	<b>215</b>
20626-250	108	108	44	74	<b>250</b>
20626-285*	113	113	46	77	<b>285</b>
20626-325	118	118	48	80	<b>325</b>
20626-365*	123	123	49	83	<b>365</b>
20626-410	128	128	51	87	<b>410</b>
20626-460*	133	133	53	90	<b>460</b>
20626-510	138	138	54	93	<b>510</b>
20626-565*	143	143	56	96	<b>565</b>

Même® TXP – Extra High Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
	A	C	B	D	
20627-095*	73	73	38	55	<b>95</b>
20627-115	78	78	41	59	<b>115</b>
20627-140*	83	83	43	62	<b>140</b>
20627-165	88	88	45	66	<b>165</b>
20627-195*	93	93	47	69	<b>195</b>
20627-225	98	98	49	73	<b>225</b>
20627-260*	103	103	51	77	<b>260</b>
20627-300	108	108	53	80	<b>300</b>
20627-345*	113	113	55	84	<b>345</b>
20627-390	118	118	57	87	<b>390</b>
20627-440*	123	123	59	91	<b>440</b>
20627-495	128	128	61	94	<b>495</b>
20627-555*	133	133	63	98	<b>555</b>
20627-615	138	138	65	101	<b>615</b>

Approximate measurements. Further sizes will be produced on demand.  
Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk \*.

## Information on Material

POLYTECH Health & Aesthetics implants are filled with a highly cross-linked gel, which stands out due to its natural consistency. Because of the "memory effect" of the gel, it always returns to its original shape after mechanical impact.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is textured. This texture is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates.

## Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

## Important Information

The implants are delivered individually packed and sterile.

**Attention:** Please note the detailed instructions for use included with each product.

POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

☎ +49 (0) 60 71 98 63 0 • 📠 +49 (0) 60 71 98 63 30

www.polytech-health-aesthetics.com • e-mail: info@polytechhealth.com

