

Optimam®

textured, silicone gel-filled
mammary implants, anatomical shape,
available in four projections

POLYTECH

Health & Aesthetics







The silicone gel-filled Optimam® implants with textured surface have a vertically oval base and an anatomical projection whose maximum point is located in the lower half of the oval base.





The consistency of these high-quality mammary implants is very close to the haptics and motional properties of the mammae. The implants are available in four projections:





Low Profile (TLP), Moderate Profile (TMP), High Profile (THP) and Extra High Profile (TXP).





The new curve measurement D, which describes the distance between the highest point of projection and the caudal implant edge, establishes a reference parameter for the correct placement of the areola in relation to the inframammary fold.

Optimam®, CE 0483

Optimam® TLP – Low Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
A	C	B	D		
20674-100*	93	107	25	43	100
20674-120	98	113	26	45	120
20674-135*	103	119	27	47	135
20674-155	108	125	28	49	155
20674-175*	113	131	28	51	175
20674-200	118	136	29	54	200
20674-220*	123	142	30	56	220
20674-250	128	148	31	58	250
20674-275*	133	154	32	60	275
20674-305	138	160	33	62	305
20674-340*	143	165	34	64	340
20674-375	148	171	34	66	375
20674-410*	153	177	35	68	410
20674-450	158	183	36	70	450

Optimam® TMP – Moderate Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
A	C	B	D		
20675-115	88	102	31	46	115
20675-135*	93	107	32	49	135
20675-155	98	113	34	51	155
20675-180*	103	119	35	53	180
20675-205	108	125	36	56	205
20675-235*	113	131	37	58	235
20675-265	118	136	39	60	265
20675-300*	123	142	40	63	300
20675-335	128	148	41	65	335
20675-370*	133	154	42	67	370
20675-415	138	160	44	70	415
20675-460*	143	165	45	72	460
20675-505	148	171	46	75	505
20675-555*	153	177	47	77	555
20675-610	158	183	49	79	610

Optimam® THP – High Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
A	C	B	D		
20676-140	88	102	38	52	140
20676-165*	93	113	39	54	165
20676-195	98	98	41	57	195
20676-225*	103	119	43	60	225
20676-255	108	125	44	62	255
20676-290*	113	131	46	65	290
20676-330	118	136	48	68	330
20676-375*	123	142	49	70	375
20676-420	128	148	51	73	420
20676-470*	133	154	53	76	470
20676-520	138	160	54	78	520
20676-580*	143	165	56	81	580

Optimam® TXP – Extra High Profile, textured					
Order no.	Dimensions (mm)				Volume (ml)
					
A	C	B	D		
20677-100*	73	84	38	48	100
20677-120	78	90	41	51	120
20677-145*	83	96	43	54	145
20677-170	88	102	45	57	170
20677-200*	93	107	47	60	200
20677-230	98	113	49	63	230
20677-265*	103	119	51	66	265
20677-305	108	125	53	69	305
20677-350*	113	131	55	72	350
20677-400	118	136	57	75	400
20677-450*	123	142	59	78	450
20677-505	128	148	61	81	505
20677-565*	133	154	63	84	565
20677-630	138	160	65	87	630

Approximate measurements. Further sizes will be produced on demand.
Silicone gel-filled, smooth sizers (non sterile, re-sterilisable) are available for products marked with an asterisk *.

Information on Material

POLYTECH Health & Aesthetics implants are filled with a highly cross-linked gel, which stands out due to its natural consistency. Because of the "memory effect" of the gel, it always returns to its original shape after mechanical impact.

The chemically and mechanically resistant, elastic shell of the implants consists of several layers of a silicone elastomer. In order to maintain a high level of safety, the shell of the silicone gel-filled implants is equipped with an additional diffusion barrier layer.

The surface of the implant shell is textured. This texture is known for its reliability and has been proven to contribute to a reduction of capsular contracture rates.

Indications

- Breast reconstruction after mastectomy or ablatio mammae
- Contour augmentation and correction
- Adjustment of asymmetries, aplasia or hypoplasia
- Combined breast and thoracic deformities
- Ptosis or atrophy
- Implant exchange

Important Information

The implants are delivered individually packed and sterile.

Attention: Please note the detailed instructions for use included with each product.

POLYTECH Health & Aesthetics GmbH

Altheimer Str. 32 • 64807 Dieburg • Germany

☎ +49 (0) 60 71 98 63 0 • 📠 +49 (0) 60 71 98 63 30

www.polytech-health-aesthetics.com • e-mail: info@polytechhealth.com

